



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® Syringe Filter, CME, non-sterile, 3.0, Ø 30 mm

€119.00
plus VAT &
Shipping

Product Images



Description

This non-sterile **qpore® syringe filter** has a hydrophilic membrane made of CME with very good thermal and chemical resistance. With an excellent flow rate, it is optimally suited for clear filtration.

- Filter housing made of polypropylene, maximum pressure resistance 5.0 bar
- Fast filtration

Additional Information

No.	6-0142
Manufacturer (Brand)	qpore
EAN	4058072173609
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	30 mm
TBST MAX	90 °C
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophilic
Material membrane	Mixed cellulose ester (MCE)
Pore size	3.0 µm
Type Connection Output	Luer lock male
Type Connection Input	Luer lock female
Area diaphragm	4.6 cm ²
Type filter	Syringe pre-filter
for medium	Liquids

