



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® syringe prefilter, CME, glass fiber prefilter, non-sterile, 0.45, Ø 30 mm

€119.00
plus VAT &
Shipping

Product Images



Description

This **non-sterile qpore® syringe filter** has a **hydrophilic membrane made of CME** with a very good thermal and chemical resistance. Due to an excellent flow rate, it is optimally suited for clear filtration. The effective filtration area of this syringe filter is 4.90 cm². The stable filter housing made of polypropylene is pressure resistant up to a maximum of 5.0 bar, which enables fast filtration.

Additional Information

No.	6-0079
Manufacturer (Brand)	qpore
EAN	4058072185626
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
Filter properties	suitable for HPLC with glass fiber prefilter
MAX operating pressure	6 bar abs.
DM outside	30 mm
Area diaphragm	4.6 cm ²
Fluid behavior	hydrophilic
Material membrane	Mixed cellulose ester (MCE)
Pore size	0.45 µm
TBST MAX	100 °C
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids

