



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## qpore® Syringe Filter, PTFE, sterile, 0.22 µm, Ø 25 mm

**€169.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The sterile hydrophobic **qpore® syringe filters** made of polytetrafluoroethylene PTFE are ideally suited for the filtration of media where high chemical resistance is required. The PTFE membranes are biologically and chemically inert and exhibit high flow rates.

- Low dead volume
- Stable at pH 1-14
- DNA-, DNase-, RNase-, pyrogen-free
- Luer connections: Luer lock female, luer lock male
- No risk of mix-up due to labeling (membrane type, pore size)
- The syringe filters are individually sterile packed to 100 pieces.

## Additional Information

No.	6-0129
Manufacturer (Brand)	qpore
EAN	4058072173111
Transport temperature	Room temperature
Color	orange
Material	Polypropylene (PP)
sterile	Yes
suitable for	Syringes
DM outside	25 mm
TBST MAX	130 °C
Filter properties	suitable for HPLC single packed
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophobic
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	3.9 cm <sup>2</sup>
Type filter	Syringe pre-filter
for medium	Liquids

