



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## qpore® Syringe Filter, PVDF, glass fiber prefilter, non-sterile, hydrophobic, 0.22, Ø 30 mm

**€129.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The non-sterile hydrophobic **qpore® syringe filters** made of polyvinylidene fluoride PVDF bind far less protein than membrane filters made of nylon, PTFE or cellulose mixed ester. They have a very high flux rate compared to other membrane materials. Due to the low protein and low drug binding properties, syringe filters made of PVDF are excellently suited for filtration in the life science sector.

## Additional Information

---

No.	6-0074
Manufacturer (Brand)	qpore
EAN	4058072173951
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	30 mm
TBST MAX	100 °C
Filter properties	suitable for HPLC
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophobic
Material membrane	Polyvinylidene fluoride (PVDF)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	4.6 cm <sup>2</sup>
Type filter	Syringe pre-filter
for medium	Liquids

