



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® Syringe Filter, PVDF, glass fiber prefilter, non-sterile, hydrophobic, 0.45, Ø 13 mm

€109.00
plus VAT &
Shipping

Product Images



Description

The non-sterile hydrophobic **qpore® syringe filters** made of polyvinylidene fluoride PVDF bind far less protein than membrane filters made of nylon, PTFE or cellulose mixed ester. They have a very high flux rate compared to other membrane materials. Due to their low protein and low drug binding properties, syringe filters made of PVDF are ideally suited for filtration in the life science sector.

- Very high flow rate
- Low protein and drug binding property
- Excellent for filtration in the life science sector

Additional Information

No.	6-0115
Manufacturer (Brand)	qpore
EAN	4058072173883
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	13 mm
TBST MAX	100 °C
Filter properties	suitable for HPLC
MAX operating pressure	6 bar abs.
Area diaphragm	1.09 cm ²
Fluid behavior	hydrophobic
Material membrane	Polyvinylidene fluoride (PVDF)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids

