



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



qpore® Syringe Filter, PVDF, sterile, hydrophilic, 0.45 µm, Ø 33 mm

€179.00
plus VAT &
Shipping

Product Images



Description

The sterile hydrophilic **qpore® syringe filters** made of polyvinylidene fluoride PVDF bind far less protein than syringe filters made of nylon, PTFE or cellulose mixed ester. They have a very high flux rate compared to other membrane materials. Due to the low protein and low drug binding property, the membrane filters made of PVDF are excellent for filtration in life science applications.

Membrane diameter 33 mm

Low dead volume

Stable at pH 1-14

DNA-, DNase-, RNase-, pyrogen-free

Luer connections: Luer lock female, Luer lock male

No risk of mix-up due to labeling (membrane type, pore size)

The syringe filters are individually sterile packed to 100 pieces.

Additional Information

No.	6-0071
Manufacturer (Brand)	qpore
EAN	4058072171773
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	Yes
suitable for	Syringes
DM outside	33 mm
TBST MAX	100 °C
Filter properties	with special chemical resistance
MAX operating pressure	6 bar abs.
Area diaphragm	4.6 cm ²
Fluid behavior	hydrophilic
Material membrane	Polyvinylidene fluoride (PVDF)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Type filter	Syringe pre-filter
for medium	Liquids

