



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## qpore® nylon syringe filter, non-sterile, 0.45 $\mu\text{m}$ , $\varnothing$ 30 mm, 100 pcs/pack

**€119.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**qpore®** offers a wide range of high quality syringe filters for various filtration applications in your laboratory. All filters are manufactured under the highest quality standards from the best raw materials.

This **non-sterile qpore® syringe filter** has a hydrophilic **nylon** membrane and is optimal for filtration of aqueous solutions and solvents. The clean and pure nylon membrane combines the fastest flow rates with low non-specific binding.

**Glass fiber pre-filter included.**

Thus, this syringe prefilter can be used wherever application or clear filtration is allowed under non-sterile conditions, as most often in sample preparation prior to HPLC or GC. The effective filtration area of the syringe filter is 1.65 or 4.90 cm<sup>2</sup>. The sturdy polypropylene filter housing is pressure resistant up to a maximum of 6.0 bar allowing fast filtration.

Features:

- Low dead volume
- Stable at pH 3-12
- Luer connections: Luer lock female, Luer cone male
- No risk of confusion due to labeling (membrane type, pore size)
- The syringe filters are packed non-sterile 100 pieces in a bag.

## Additional Information

No.	6-0019
Manufacturer (Brand)	qpore
EAN	4058072172657
Transport temperature	Room temperature
Color	orange
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	30 mm
TBST MAX	100 °C
Filter properties	suitable for HPLC
MAX operating pressure	7.5 bar abs.
Fluid behavior	hydrophilic
Material housing	Polypropylene (PP)
Material membrane	Polyamide 6.6 (PA 6.6. Nylon®)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	4.6 cm <sup>2</sup>
Type filter	Syringe pre-filter
for medium	Liquids

