



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**qpore® syringe filter, PTFE, glass fiber prefilter,  
non-sterile, hydrophobic, 0.45, Ø 13 mm**

**€109.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

**qpore®** offers a wide range of high quality syringe filters for various filtration applications in your laboratory. All filters are manufactured under the highest quality standards from the best raw materials.

This **non-sterile qpore® syringe filter** has a membrane made of **PTFE** with a high resistance to solvents, acids and bases. Thus, this syringe prefilter can be used wherever application or clear filtration is allowed under non-sterile conditions, as mostly in sample preparation prior to HPLC or GC. **Including pre-filter made of glass fiber.** The sturdy filter housing made of polypropylene is pressure resistant up to a maximum of 6.0 bar which allows fast filtration.

Features:

- Low dead volume
- Stable at pH 1-14
- Luer connections: Luer lock female, Luer cone male
- No risk of mix-up due to labeling (membrane type, pore size)
- The syringe filters are packed non-sterile 100 pieces in a bag

## Additional Information

No.	6-0113
Manufacturer (Brand)	qpore
EAN	4058072173876
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	13 mm
TBST MAX	100 °C
Filter properties	suitable for HPLC hydrophobic
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophobic
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.45 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	1.09 cm <sup>2</sup>
Type filter	Syringe pre-filter
for medium	Liquids

