



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**qpore® Syringe Filter, PTFE, glass fiber prefilter,
non-sterile, hydrophilic, 0.22 µm, Ø 25 mm**

**€119.00
plus VAT &
Shipping**

Product Images



Description

The non-sterile hydrophilic **qpore® syringe filters** made of polytetrafluoroethylene PTFE are ideally suited for the filtration of media where high chemical resistance is required. The PTFE membranes are biologically and chemically inert and exhibit high flow rates.

- Low dead volume
- Stable at pH 1-14
- DNA-, DNase-, RNase-, pyrogen-free
- Luer connections: Luer lock female, luer lock male
- No risk of mix-up due to labeling (membrane type, pore size)
- The syringe filters are individually sterile packed to 100 pieces.

Additional Information

No.	6-0134
Manufacturer (Brand)	qpore
EAN	4058072177089
Transport temperature	Room temperature
Material	Polypropylene (PP)
sterile	No
suitable for	Syringes
DM outside	25 mm
TBST MAX	100 °C
Filter properties	suitable for HPLC
MAX operating pressure	6 bar abs.
Fluid behavior	hydrophilic
Material membrane	Polytetrafluoroethylene (PTFE)
Pore size	0.22 µm
Type Connection Output	Luer cone male
Type Connection Input	Luer lock female
Area diaphragm	3.9 cm ²
Type filter	Syringe pre-filter
for medium	Liquids

