



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Sartorius Vivaspin 2, MWCO 10 kDa, polyethersulfone, 25 pcs

**€174.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Vivaspin® 2** combines the speed of the classic Vivaspin® products with low internal surface area and membrane area to achieve excellent recovery from highly diluted solutions. Recovery of concentrate possible by reverse centrifugation.

**For 0.4 to 2 ml sample solution.**

Tubes fit all standard fixed rotors for 1.5 and 2.0 ml reaction tubes.

## Additional Information

No.	ST-2632
Manufacturer (Brand)	Sartorius
old neoLab article no.	8-6023
EAN	4058072370558
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	2 ml
sterile	No
Type closure/ lid	Screw cap
TBST MAX	25 °C
TBST MIN	4 °C
Filter properties	with special chemical resistance
Floor outside	conical
Color graduation	black
Format/ Packing	in rack
MWCO	10 kDa
Material membrane	Polyethersulfone (PES)
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type filter	Vertical twin diaphragm
Type Tube/ Micro Vessel	Ultrafiltration unit

