



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius Vivaspin 2, MWCO 50 kDa, polyethersulfone, 25 pcs

€174.00
plus VAT &
Shipping

Product Images



Description

With Vivaspin ultracentrifugal tubes, save your protein. With the polyethersulfone (PES) membrane, your wic

Vivaspin® 2 combines the speed of classic Vivaspin® products with low internal surface area and membrane area to achieve excellent recovery from highly diluted solutions. Recovery of concentrate possible by reverse centrifugation.

For 0.4 to 2 ml sample solution.

Tubes fit all standard fixed rotors for 1.5 and 2.0 ml reaction tubes.

Additional Information

No.	ST-2638
Manufacturer (Brand)	Sartorius
old neolab article no.	8-6027
EAN	4058072370619
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	2 ml
sterile	No
Type closure/ lid	Screw cap
TBST MAX	25 °C
TBST MIN	4 °C
Filter properties	with special chemical resistance
Floor outside	conical
Color graduation	black
MWCO	50 kDa
Material membrane	Polyethersulfone (PES)
Format/ Packing	in rack
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type filter	Vertical twin diaphragm
Type Tube/ Micro Vessel	Ultrafiltration unit

