



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius Vivaspin 2, MWCO 1000 kDa, polyethersulfone, 25 pcs

€174.00
plus VAT &
Shipping

Product Images



Description

Vivaspin® 2 combines the speed of the classic Vivaspin® products with low internal surface area and membrane area to achieve excellent recovery from highly diluted solutions. Recovery of concentrate possible by reverse centrifugation.

For 0.4 to 2 ml sample solution.

Tubes fit all standard fixed rotors for 1.5 and 2.0 ml reaction tubes.

Additional Information

| | |
|-------------------------|----------------------------------|
| No. | ST-2644 |
| Manufacturer (Brand) | Sartorius |
| old neoLab article no. | 8-6033 |
| EAN | 4058072370671 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | transparent (crystal clear) |
| Volume | 2 ml |
| sterile | No |
| Type closure/ lid | Screw cap |
| TBST MAX | 25 °C |
| TBST MIN | 4 °C |
| Filter properties | with special chemical resistance |
| Floor outside | conical |
| Color graduation | black |
| MWCO | 1000 kDa |
| Material membrane | Polyethersulfone (PES) |
| Format/ Packing | in rack |
| Material tube | Polycarbonate (PC) |
| Material closure/ lid | Polypropylene (PP) |
| Type filter | Vertical twin diaphragm |
| Type Tube/ Micro Vessel | Ultrafiltration unit |

