



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Sartorius Vivaspin 20, 0.2 $\mu\text{m}$ , polyethersulfone, 12 pcs

**€177.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

**Vivaspin® 20 centrifugal concentrators** are designed to offer increased volume flexibility and performance. Vivaspin® 20 can process samples up to 20 ml in swing-out rotors or 14 ml in fixed-angle rotors for 50 ml centrifuge tubes. The remaining volume can be easily read from the printed scale and the modified dead-stop reservoir facilitates direct pipetting of the concentrate.

**For 5 to 20 ml sample solution.**

Tubes fit all standard fixed rotors for 50 ml reaction tubes.

## Additional Information

No.	ST-2731
Manufacturer (Brand)	Sartorius
old neoLab article no.	8-6119
EAN	4066292265059
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	20 ml
sterile	No
Type closure/ lid	Screw cap
TBST MAX	25 °C
TBST MIN	4 °C
Filter properties	with special chemical resistance
Floor outside	conical
Color graduation	black
Material membrane	Polyethersulfone (PES)
Pore size	2 µm
Format/ Packing	in rack
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type filter	Vertical twin diaphragm
Type Tube/ Micro Vessel	Ultrafiltration unit

