



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Sartorius Vivaspin 6, MWCO 1000 kDa, polyethersulfone, 25 pcs

**€250.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Vivaspin® 6 centrifugal concentrators** are designed to provide increased volume flexibility and performance. Vivaspin® 6 can process up to 6 ml of sample in fixed-angle or swing-out rotors for 15 ml standard pointed-bottom tubes. The units are equipped with vertically mounted twin membranes that provide unprecedented filtration speeds and concentration factors over 100-fold. The remaining volume can be easily read from the printed scale and the modified dead-stop reservoir facilitates direct withdrawal of the concentrate with the pipette.

**For 2 to 6 ml sample solution.**

Tubes fit all common fixed rotors for 15 ml reaction tubes.

## Additional Information

No.	ST-2689
Manufacturer (Brand)	Sartorius
old neoLab article no.	8-6078
EAN	4058072370985
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	6 ml
sterile	No
Type closure/ lid	Screw cap
TBST MAX	25 °C
TBST MIN	4 °C
Filter properties	with special chemical resistance
Floor outside	conical
Color graduation	black
Format/ Packing	in rack
MWCO	1000 kDa
Material membrane	Polyethersulfone (PES)
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type filter	Vertical twin diaphragm
Type Tube/ Micro Vessel	Ultrafiltration unit

