



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius Vivaspin 2, starter pack (5 pieces each: 5 k, 10 k, 30 k, 50 k, 100 kDa) PU: 25 pcs

**€183.00
plus VAT &
Shipping**

Product Images



Description

Vivaspin® 2 combines the speed of the classic Vivaspin® products with low internal surface area and membrane area to achieve excellent recovery from highly diluted solutions. Recovery of concentrate possible by reverse centrifugation.

For 0.4 to 2 ml sample solution.

Tubes fit all standard fixed rotors for 1.5 and 2.0 ml reaction tubes.

Additional Information

No.	ST-2658
Manufacturer (Brand)	Sartorius
old neoLab article no.	8-6037
EAN	4058072370732
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	2 ml
sterile	No
Type closure/ lid	Screw cap
TBST MAX	25 °C
TBST MIN	4 °C
Filter properties	with special chemical resistance
Floor outside	conical
Color graduation	black
MWCO	5 kDa 10 kDa 30 kDa 50 kDa 100 kDa
Number of parts in set	25
Material membrane	Polyethersulfone (PES)
Format/ Packing	in rack
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type filter	Vertical twin diaphragm
Type Tube/ Micro Vessel	Ultrafiltration unit

