



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius Vivaspin 20, starter pack (5 k, 10 k, 30 k, 50 k, 100 kDa), PU: 12 pcs

**€177.00
plus VAT &
Shipping**

Product Images



Description

Vivaspin® 20 centrifugal concentrators are designed to offer increased volume flexibility and performance. Vivaspin® 20 can process samples up to 20 ml in swing-out rotors or 14 ml in fixed-angle rotors for 50 ml centrifuge tubes. The remaining volume can be easily read from the printed scale and the modified dead-stop reservoir facilitates direct pipetting of the concentrate.

For 5 to 20 ml sample solution.

Tubes fit all standard fixed rotors for 50 ml reaction tubes.

Additional Information

No.	ST-2735
Manufacturer (Brand)	Sartorius
old neoLab article no.	8-6121
EAN	4066292265066
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	20 ml
sterile	No
Type closure/ lid	Screw cap
TBST MAX	25 °C
TBST MIN	4 °C
Filter properties	with special chemical resistance
Number of parts in set	12
Floor outside	conical
Color graduation	black
Format/ Packing	in rack
MWCO	5 kDa 10 kDa 30 kDa 50 kDa 100 kDa
Material membrane	Polyethersulfone (PES)
Material tube	Polycarbonate (PC)
Material closure/ lid	Polypropylene (PP)
Type filter	Vertical twin diaphragm
Type Tube/ Micro Vessel	Ultrafiltration unit

