



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Sieve cloth Polyester Monolen, mesh size 1000 µm, 100 x 102 cm

€100.00
plus VAT &
Shipping

Product Images



Description

Polyester monolayer mesh for filtration or chromatography columns.
Operating temperature dry from -75°C to +150°C.

Additional Information

No.	2-4054
Manufacturer (Brand)	neoLab
EAN	4058072001926
Transport temperature	Room temperature
Wide	1020 mm
Length	1000 mm
TBST MAX	150 °C
TBST MIN	-75 °C
Mesh size	1000 µm
Material screen cloth/insert	Polyester (PEs)
Format/ Packing	Roll

