



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® sieve cloth made of polyester monols,  
mesh size 500 µm, 115 x 100 cm**

**€53.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

**Polyester monolayer mesh** for filtration or chromatography columns.  
Operating temperature dry from -75°C to +150°C.

## Additional Information

---

No.	2-4059
Manufacturer (Brand)	neolab
EAN	4058072083427
Transport temperature	Room temperature
Wide	1150 mm
Length	1000 mm
TBST MAX	150 °C
TBST MIN	-75 °C
Format/ Packing	Roll
Mesh size	500 µm
Material screen cloth/insert	Polyester (PEs)

