



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Sieve cloth Polyester Monolen, mesh size 400 µm, 100 x 102 cm

**€55.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Polyester monolayer mesh** for filtration or chromatography columns.  
Operating temperature dry from -75°C to +150°C.

## Additional Information

---

No.	2-4060
Manufacturer (Brand)	neoLab
EAN	4058072001933
Transport temperature	Room temperature
Wide	1020 mm
Length	1000 mm
Format/ Packing	Roll
Mesh size	400 µm
Material screen cloth/insert	Polyester (PEs)
TBST MAX	150 °C
TBST MIN	-75 °C

