



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Sieve cloth Polyester Monolen, mesh size 200 µm, 100 x 110 cm

**€51.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Polyester monolayer mesh** for filtration or chromatography columns.  
Operating temperature dry from -75°C to +150°C.

## Additional Information

---

No.	2-4063
Manufacturer (Brand)	neoLab
EAN	4058072001964
Transport temperature	Room temperature
Wide	1020 mm
Length	1000 mm
TBST MAX	150 °C
TBST MIN	-75 °C
Format/ Packing	Roll
Mesh size	200 µm
Material screen cloth/insert	Polyester (PEs)

