



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Sieve cloth Polyester Monolen, mesh size 40 µm, 100 x 102 cm

€168.20
plus VAT &
Shipping

Product Images



Description

Polyester monolayer mesh for filtration or chromatography columns.
Operating temperature dry from -75°C to +150°C.

Additional Information

| | |
|------------------------------|------------------|
| No. | 2-4071 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072083441 |
| Transport temperature | Room temperature |
| Wide | 1020 mm |
| Length | 1000 mm |
| TBST MAX | 150 °C |
| TBST MIN | -75 °C |
| Format/ Packing | Roll |
| Mesh size | 40 µm |
| Material screen cloth/insert | Polyester (PEs) |

