





neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® sieve cloth PP monofil, mesh size 3360
µm, 100 x 102 cm**

**€93.60
plus VAT &
Shipping**



Description

Monofilament polypropylenefabric with excellent chemical resistance.

For demanding technical applications.

Operating temperature -30°C to +90°C dry.

Moisture absorption < 0.01%.

Additional Information

No.	4-1430
Manufacturer (Brand)	neoLab
EAN	4058072022266
Transport temperature	Room temperature
Wide	1020 mm
Length	1000 mm
Format/ Packing	Roll
Mesh size	3360 µm
Material screen cloth/insert	Polypropylene (PP)
TBST MAX	90 °C
TBST MIN	-30 °C

