




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® sieve cloth made of PP monofilament,
mesh size 1000 µm, 100 x 102 cm**

**€77.30
plus VAT &
Shipping**



Description

Monofilament polypropylene fabric with excellent chemical resistance.
For demanding technical applications.
Operating temperature -30°C to +90°C dry.
Moisture absorption < 0.01%.

Additional Information

No.	4-1435
Manufacturer (Brand)	neoLab
EAN	4058072083984
Transport temperature	Room temperature
Wide	1020 mm
Length	1000 mm
TBST MAX	90 °C
TBST MIN	-30 °C
Format/ Packing	Roll
Mesh size	1000 µm
Material screen cloth/insert	Polypropylene (PP)

