



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® sieve cloth PP monofil, mesh size 105
µm, 100 x 102 cm**

**€56.90
plus VAT &
Shipping**

Product Images



Description

Monofilament polypropylene fabric with excellent chemical resistance.
For demanding technical applications.
Operating temperature -30°C to +90°C dry.
Moisture absorption < 0.01%.

Additional Information

No.	4-1444
Manufacturer (Brand)	neoLab
EAN	4058072084042
Transport temperature	Room temperature
Wide	1020 mm
Length	1000 mm
TBST MAX	90 °C
TBST MIN	-30 °C
Mesh size	105 µm
Material screen cloth/insert	Polypropylene (PP)
Format/ Packing	Roll

