



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® analytical funnel, d = 65 mm

€36.66 As low
as
€28.20
plus VAT &
Shipping

Product Images



Description

DURAN® Analytical funnel, Ø 65 mm Transferring liquids and solids is a common working method in which the funnel is used. Due to its shape, the funnel can hold a relatively large quantity of liquid or solid and transfer it through a constriction into a vessel with a smaller(er) opening. Funnels come in different sizes and designs, which are selected according to their field of application. Special features and advantages: For rapid filtration Ideally suited for rapid filtration using paper round filters due to inserted glass webs Burned-in and therefore permanently durable white printed image Application: Transferring and filtering substances, Spilling of substances is avoided, Filtration

Additional Information

No.	SD-0480
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112133137
EAN	4058072177478
Manufacturer EAN	4032051008663
Gross weight	0.088
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.06300 kg
diameter	65 mm
diameter stem outside	9 mm
Material	Borosilicate glass 3.3
Length of stem	150 mm
Country of origin	Germany
sterile	No
Region of origin	Baden-Württemberg
Packaging volume	0.00827800 cc
Customs tariff number	70172000

