



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® analytical funnel, d = 80 mm

€43.21 As low as
€33.24
plus VAT &
Shipping

Product Images



Description

DURAN® Analytical funnel, Ø 80 mm Transferring liquids and solids is a common working method in which the funnel is used. Due to its shape, the funnel can hold a relatively large quantity of liquid or solid and transfer it through a constriction into a vessel with a smaller(er) opening. Funnels come in different sizes and designs, which are selected according to their field of application. Special features and advantages: For rapid filtration Ideally suited for rapid filtration using paper round filters due to inserted glass webs Burned-in and therefore permanently durable white printed image Application: Transferring and filtering substances, Spilling of substances is avoided, Filtration

Additional Information

No.	SD-0481
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112133141
EAN	4058072188405
Manufacturer EAN	4032051008670
Gross weight	0.137
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.10900 kg
diameter	80 mm
diameter stem outside	9 mm
Material	Borosilicate glass 3.3
Length of stem	150 mm
Country of origin	Germany
sterile	No
Region of origin	Baden-Württemberg
Packaging volume	0.01090700 cc
Customs tariff number	70172000

