



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Powder funnel, short and wide handle, d = 80 mm

€22.48 As low
as
€17.29
plus VAT &
Shipping

Product Images



Description

DURAN® Powder funnel, Ø 80 mm Transferring liquids and solids is a common working method in which the funnel is used. Due to its shape, the funnel can hold a relatively large quantity of liquid or solid and transfer it through a constriction into a vessel with a smaller(er) opening. Funnels come in different sizes and designs, which are selected according to their field of application. Special features and advantages: With short and wide stem Suitable for transferring solids due to large stem diameter Burnt and therefore durable white printed image Application: Transferring and filtering substances, Spillage of substances is avoided, Transfer

Additional Information

No.	SD-0514
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112135441
EAN	4058072182885
Manufacturer EAN	4032051009103
Gross weight	0.073
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.05400 kg
diameter	80 mm
diameter stem outside	24 mm
Material	Borosilicate glass 3.3
Length of stem	30 mm
Country of origin	Hungary
sterile	No
Packaging volume	0.00799300 cc
Customs tariff number	70172000

