



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz filter push, inner Ø 50 mm, height 132 mm

€25.57
plus VAT &
Shipping

Product Images



Description

Filter lugs made of borosilicate glass 3.3
Suitable for holding filter crucibles

Additional Information

No.	LZ-2283
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4287
EAN	4058072320164
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	50 mm
Color	transparent (clear)
Material	Borosilicate glass
sterile	No
Form	cylindrical
Height	132 mm

