



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® glass funnel 45 mm Ø, stem length 45 mm

€4.46
plus VAT &
Shipping

Product Images



Description

Angle 60°, with short handle.

Additional Information

No.	E-7043
Manufacturer (Brand)	neoLab
EAN	4058072014841
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	45 mm
diameter stem outside	3 mm
Material	Lime soda glass
Length of stem	45 mm
sterile	No
Height	80 mm
TBST MAX	450 °C
Funnel features	with special chemical resistance
DM neck inside	3 mm
Length neck	45 mm

