



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




## Hecht glass funnel 30 mm Ø, stem length 30 mm

**€2.24** As low as

**€1.96**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Angle 60°, with short handle.

## Additional Information

---

No.	WQ-2884
Manufacturer (Brand)	Hecht
old neoLab article no.	E-7040
EAN	4058072443849
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	30 mm
diameter stem outside	6 mm
Material	Lime soda glass
Length of stem	30 mm
sterile	No
Height	40 mm
TBST MAX	450 °C
Funnel features	with special chemical resistance
DM neck inside	6 mm
Length neck	30 mm

