



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht glass funnel 40 mm Ø, stem length 40 mm

€2.73 As low as
€2.34
plus VAT &
Shipping

Product Images



Description

Angle 60°, with short handle.

Additional Information

No.	WQ-2886
Manufacturer (Brand)	Hecht
old neoLab article no.	E-7042
EAN	4058072443856
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	40 mm
diameter stem outside	8 mm
Material	Lime soda glass
Length of stem	40 mm
sterile	No
Height	50 mm
TBST MAX	450 °C
Funnel features	with special chemical resistance
DM neck inside	8 mm
Length neck	40 mm

