



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Urbanti funnel 100 mm Ø, 7 mm Ø drain, 4 pcs./bag

€38.60
plus VAT &
Shipping

Product Images



Description

Unsurpassed for rapid filtration under non-vacuum.

With spiral ribs on the inside to increase the permeability of filter paper.

Vertical ribs on the outer wall to prevent pressure build-up.

Autoclavable, heat resistant up to 177°C, chemically resistant to corrosive liquids.

Additional Information

No.	6-1113
Manufacturer (Brand)	neolab
EAN	4058072024789
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	111 mm
diameter stem outside	10 mm
Material	Polymethylpenten (PMP, TPX®)
Length of stem	108 mm
sterile	No
Height	194 mm
TBST MAX	177 °C
Funnel features	with special chemical resistance inside with spiral ribs
for filters with DM	150 mm 185 mm
DM neck inside	7 mm
Length neck	108 mm

