



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Urbanti funnel 196 mm Ø, 20 mm Ø drain, 2 pcs./bag

€110.00
plus VAT &
Shipping

Product Images



Description

Unsurpassed for rapid filtration under non-vacuum.

With spiral ribs on the inside to increase the permeability of filter paper.

Vertical ribs on the outer wall to prevent pressure build-up.

Autoclavable, heat resistant up to 177°C, chemically resistant to corrosive liquids.

Additional Information

No.	6-1115
Manufacturer (Brand)	neolab
EAN	4058072024802
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	209 mm
diameter stem outside	26 mm
Material	Polymethylpenten (PMP, TPX®)
Length of stem	160 mm
sterile	No
Height	316 mm
TBST MAX	177 °C
Funnel features	with special chemical resistance inside with spiral ribs
for filters with DM	320 mm 385 mm
DM neck inside	20 mm
Length neck	160 mm

