



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnels acc. to Gilson, with PTFE
plug, 2,000 ml, NS 29/32, bore 6.0 mm**

**€220.90
plus VAT &
Shipping**

Product Images



Description

With NS sleeve and plastic plug (PE), without pitch.

With PTFE plug, with threaded stopcock safety device, with labeling field.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0859
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2859
EAN	4058072305994
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	2000 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Material closure/ lid	Polyethylene (PE)

