



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz separatory funnel acc. to Squibb, with solid glass plug, 50 ml, NS 19/26, bore 2.5 mm, graduation 2 ml**

**€72.92** As low as  
**€56.09**  
**plus VAT & Shipping**

## Product Images



## Description

---

**Separatory funnel** acc. to Squibb, with glass plug, with graduation.

With NS sleeve and plastic stopper (PE), with solid NS glass plug, with threaded stopcock safety device.

With graduation for reading the Zirka volume.

Material: borosilicate glass 3.3.

## Additional Information

---

No.	LZ-0808
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2808
EAN	4058072305482
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Borosilicate glass
Standard ground	19/26
Division	2 ml
Material closure/ lid	Polyethylene (PE)

