



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz separatory funnel acc. to Squibb, with solid glass plug, 100 ml, NS 19/26, bore 2.5 mm, graduation 5 ml

€68.28 *As low*
as
€52.52
**plus VAT &
Shipping**

Product Images



Description

Separatory funnel acc. to Squibb, with glass plug, with graduation.

With NS sleeve and plastic stopper (PE), with solid NS glass plug, with threaded stopcock safety device.

With graduation for reading the Zirka volume.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0809
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2809
EAN	4058072305499
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Borosilicate glass
Standard ground	19/26
Division	5 ml
Material closure/ lid	Polyethylene (PE)

