



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz separatory funnel acc. to Squibb, with solid glass plug, 500 ml, NS 29/32, bore 4.0 mm, without graduation

€84.85 As low
as
€65.27
plus VAT &
Shipping



Description

Separatory funnel acc. to Squibb, with glass plug, without graduation.

With NS sleeve and plastic stopper (PE).

With solid NS glass plug, with threaded stopcock.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0805
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2805
EAN	4058072305451
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Borosilicate glass
Material closure/ lid	Polyethylene (PE)
Standard ground	29/32

