



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz separatory funnel acc. to Squibb, with PTFE plug, 500 ml, NS 29/32, bore 4.0 mm, graduation 10 ml**

**€105.11** As low  
as  
**€80.85**  
**plus VAT &  
Shipping**



## Description

---

**Separatory funnels** acc. to Squibb, with PTFE plug, with graduation.

With NS sleeve and plastic stopper (PE), with threaded stopcock.

With graduation for reading off the Zirka volume.

Material: borosilicate glass 3.3.

## Additional Information

---

No.	LZ-0823
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2823
EAN	4058072305635
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Division	10 ml
Material closure/ lid	Polyethylene (PE)

