



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnels acc. to Squibb, with
PTFE plug, 500 ml, NS 29/32, bore 4.0 mm,
without graduation**

€97.62 As low
as
€75.09
plus VAT &
Shipping



Description

Separatory funnels acc. to Squibb, with PTFE plug, without graduation.
With NS sleeve and plastic stopper (PE), with threaded stopcock.
Material: borosilicate glass 3.3

Additional Information

No.	LZ-0817
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2817
EAN	4058072305574
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Material closure/ lid	Polyethylene (PE)

