



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnels acc. to Squibb, with
PTFE plug, 100 ml, NS 19/26, bore 2.5 mm, amber
glass, without graduation**

€106.46 As low
as
€81.89
plus VAT &
Shipping



Description

Separatory funnels acc. to Squibb, amber glass, without graduation, with PTFE plug
With NS sleeve and plastic stopper (PE), with threaded stopcock.
Material: borosilicate glass 3.3

Additional Information

No.	LZ-0833
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2833
EAN	4058072305734
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Borosilicate glass
Volume	100 ml
sterile	No
Form	pear-shaped
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	19/26
Material closure/ lid	Polyethylene (PE)

