



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz separatory funnels acc. to Squibb, with PTFE plug, 250 ml, NS 29/32, bore 4.0 mm, amber glass, without graduation**

**€127.54** As low  
as  
**€98.11**  
**plus VAT &  
Shipping**



## Description

---

Separatory funnels acc. to Squibb, amber glass, without graduation, with PTFE plug  
With NS sleeve and plastic stopper (PE), with threaded stopcock.  
Material: borosilicate glass 3.3

## Additional Information

---

|                        |  |
|------------------------|--|
| No.                    | LZ-0834  |
| Manufacturer (Brand)   | Lenz   |
| old neoLab article no. | B-2834   |
| EAN                    | 4058072305741  |
| Transport temperature  | Room temperature   |
| autoclavable at 121°C  | Yes  |
| Color                  | brown  |
| Material               | Borosilicate glass   |
| Volume                 | 250 ml   |
| sterile                | No   |
| Form                   | pear-shaped  |
| Type closure/ lid      | NS plug  |
| Funnel features        | with closure/lid with standard ground joint with NS stopcock/valve |
| Material chucks        | Polytetrafluoroethylene (PTFE)                                     |
| Standard ground        | 29/32  |
| Material closure/ lid  | Polyethylene (PE)  |

