



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnel, conical, with solid glass
plug, 500 ml, NS 29/32, bore 4.0 mm, graduation
10 ml**

**€146.04
plus VAT &
Shipping**

Product Images



Description

Conical **separating funnel**, with glass plug, with graduation.

With NS sleeve and plastic stopper (PE), with graduation for reading the Zirka volume.

With solid NS glass plug, with thread lock.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0775
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2775
EAN	4058072305154
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Borosilicate glass
Standard ground	29/32
Division	10 ml
Material closure/ lid	Polyethylene (PE)

