



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz separatory funnel, conical, with solid glass  
plug, 1,000 ml, NS 29/32, bore 6.0 mm,  
graduation 20 ml**

**€210.20  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Conical **separating funnel**, with glass plug, with graduation.

With NS sleeve and plastic stopper (PE), with graduation for reading the Zirka volume.

With solid NS glass plug, with thread lock.

Material: borosilicate glass 3.3.

## Additional Information

---

No.	LZ-0776
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2776
EAN	4058072305161
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Borosilicate glass
Standard ground	29/32
Division	20 ml
Material closure/ lid	Polyethylene (PE)

