



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Separatory Funnel, conical, with solid glass
plug, 100 ml, NS 19/26, bore 2.5 mm, without
graduation**

**€101.00
plus VAT &
Shipping**

Product Images



Description

Conical separating funnel, with glass plug without graduation.

With NS sleeve and plastic stopper (PE), without graduation, with solid NS glass plug, with threaded stopcock safety device. Made of borosilicate glass 3.3.

Additional Information

No.	LZ-0767
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2767
EAN	4058072305079
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Borosilicate glass
Standard ground	19/26
Material closure/ lid	Polyethylene (PE)

