



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnel, conical, with solid glass
plug, 500ml, NS 29/32, bore 4.0 mm, amber
glass, without graduation**

**€232.84
plus VAT &
Shipping**

Product Images



Description

Conical **separating funnel**, amber glass, without graduation, with glass plug.
With NS sleeve and plastic stopper (PE), without graduation, with solid NS glass plug, with thread lock.
Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0793
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2793
EAN	4058072305338
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Borosilicate glass
Volume	500 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Borosilicate glass
Standard ground	29/32
Material closure/ lid	Polyethylene (PE)

