



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz separatory funnel, conical, with PTFE plug, 50 ml, NS 19/26, bore 2.5 mm, graduation 2 ml

€112.14
plus VAT &
Shipping

Product Images



Description

Separating funnel, conical, with PTFE plug, with graduation.

With NS sleeve and plastic stopper (PE), with graduation for reading off the Zirka volume.

With PTFE plug, with threaded stopcock.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0784
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2784
EAN	4058072305246
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Polytetrafluoroethylene (PTFE)
Material closure/ lid	Polyethylene (PE)
Standard ground	19/26
Division	2 ml

