



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnels, conical, with PTFE plug,
1,000 ml, NS 29/32, bore 6.0 mm, without
graduation**

**€215.60
plus VAT &
Shipping**

Product Images



Description

Separatory funnel, conical, with PTFE plug, without graduation.

With NS sleeve and plastic stopper (PE), without graduation, with PTFE plug, with threaded stopcock.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0782
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2782
EAN	4058072305222
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Material closure/ lid	Polyethylene (PE)

