



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz Separatory Funnel, conical, with PTFE plug,  
500 ml, NS 29/32, bore 4.0 mm, amber glass,  
without graduation**

**€251.28  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Conical **separating funnel**, amber glass, without graduation, with PTFE plug.  
With NS sleeve and plastic stopper (PE), with threaded stopcock.  
Material: borosilicate glass 3.3.

## Additional Information

---

No.	LZ-0799
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2799
EAN	4058072305390
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Borosilicate glass
Volume	500 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Material closure/ lid	Polyethylene (PE)

