



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Separatory Funnel, conical, with PTFE plug,
1.000 ml, NS 29/32, bore 6.0 mm, amber glass,
without graduation**

**€326.96
plus VAT &
Shipping**

Product Images



Description

Conical **separating funnel**, amber glass, without graduation, with PTFE plug.
With NS sleeve and plastic stopper (PE), with threaded stopcock.
Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0800
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2800
EAN	4058072305406
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Form	conical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Material closure/ lid	Polyethylene (PE)

