



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz separatory funnel, cylindrical, with solid glass plug, 50 ml, NS 29/32, bore 2.5 mm, graduation 2 ml

**€108.54
plus VAT &
Shipping**

Product Images



Description

Cylindrical shape with smooth drain tube, with NS sleeve and plastic stopper (PE).

With solid NS glass plug, with threaded stopcock safety device, with graduation for reading off the Zirka volume.

Material: Borosilicate 3.3

Additional Information

No.	LZ-0838
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2838
EAN	4058072305789
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
sterile	No
Form	cylindrical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Borosilicate glass
Standard ground	29/32
Division	2 ml
Material closure/ lid	Polyethylene (PE)

